

PTP IA-2β rabbit pAb

Cat No.: ES6943

For research use only

Overview

Product Name PTP IA-2β rabbit pAb

Host species Rabbit

Applications WB;ELISA;IHC Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000 **Immunogen** The antiserum was produced against synthesized

peptide derived from human PTPRN2. AA

range:206-255

Specificity PTP IA-2β Polyclonal Antibody detects endogenous

levels of PTP IA-2β protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Receptor-type tyrosine-protein phosphatase N2

Gene Name PTPRN2

Cellular localization Cytoplasmic vesicle, secretory vesicle membrane;

Single-pass type I membrane protein . Cytoplasmic

vesicle, secretory vesicle, synaptic vesicle

membrane; Single-pass type I membrane protein. Predominantly found on dense-core secretory

granules. Sorting

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 111kD
Human Gene ID 5799
Human Swiss-Prot Number Q92932

Alternative Names PTPRN2; KIAA0387; Receptor-type tyrosine-protein

phosphatase N2; R-PTP-N2; Islet cell

autoantigen-related protein; IAR; ICAAR; Phogrin


Background This gene encodes a protein with sequence




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com



similarity to receptor-like protein tyrosine phosphatases. However, tyrosine phosphatase activity has not been experimentally validated for this protein. Studies of the rat ortholog suggest that the encoded protein may instead function as a phosphatidylinositol phosphatase with the ability to dephosphorylate phosphatidylinositol 3-phosphate and phosphatidylinositol 4,5-diphosphate, and this function may be involved in the regulation of insulin secretion. This protein has been identified as an autoantigen in insulin-dependent diabetes mellitus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2015],



Western blot analysis of lysates from NIH/3T3 cells, using PTPRN2 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human Breast cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

