

Peroxin 19 rabbit pAb

Cat No.: ES6950

For research use only

Overview

Product Name Peroxin 19 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human PEX19. AA

range:219-268

Specificity Peroxin 19 Polyclonal Antibody detects endogenous

levels of Peroxin 19 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Peroxisomal biogenesis factor 19

Gene Name PEX19

Cytoplasm . Peroxisome membrane ; Lipid-anchor ;

Cytoplasmic side . Mainly cytoplasmic. Some fraction membrane-associated to the outer surface

of peroxisomes...

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 33kD
Human Gene ID 5824
Human Swiss-Prot Number P40855

Alternative Names PEX19; HK33; PXF; OK/SW-cl.22; Peroxisomal

biogenesis factor 19; 33 kDa housekeeping protein;

Peroxin-19; Peroxisomal farnesylated protein peroxisomal biogenesis factor 19(PEX19) Homo




Background


+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com



sapiens This gene is necessary for early peroxisomal biogenesis. It acts both as a cytosolic chaperone and as an import receptor for peroxisomal membrane proteins (PMPs). Peroxins (PEXs) are proteins that are essential for the assembly of functional peroxisomes. The peroxisome biogenesis disorders (PBDs) are a group of genetically heterogeneous autosomal recessive, lethal diseases characterized by multiple defects in peroxisome function. These disorders have at least 14 complementation groups, with more than one phenotype being observed for some complementation groups. Although the clinical features of PBD patients vary, cells from all PBD patients exhibit a defect in the import of one or more classes of peroxisomal matrix proteins into the organelle. Defects in this gene are a cause of Zellweger syndrome (ZWS), as well as peroxisome biogenesis disorder complementation group 14 (PBD-CG14), which is



Western blot analysis of lysates from MCF-7 cells, using PEX19 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using PEX19 antibody.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com





Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).



+86-27-59760950