

## S-100A3 rabbit pAb

## Cat No.:ES7113

For research use only

## Overview

Product Name	S-100A3 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/40000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
_	peptide derived from human S100A3. AA
	range:26-75
Specificity	S-100A3 Polyclonal Antibody detects endogenous
	levels of S-100A3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Protein S100-A3
Gene Name	S100A3
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	22kD
Human Gene ID	6274
Human Swiss-Prot Number	P33764
Alternative Names	S100A3; S100E; Protein S100-A3; Protein S-100E;
	S100 calcium-binding protein A3
Background	S100 calcium binding protein A3(S100A3) Homo
	sapiens The protein encoded by this gene is a
	member of the S100 family of proteins containing 2
	EF-hand calcium-binding motifs. S100 proteins are
	localized in the cytoplasm and/or nucleus of a wide
	range of cells, and involved in the regulation of a



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein has the highest content of cysteines of all S100 proteins, has a high affinity for Zinc, and is highly expressed in human hair cuticle. The precise function of this protein is unknown. [provided by RefSeq, Jul 2008],

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using S100A3 Antibody. The picture on the right is blocked with the synthesized peptide.





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C