

Sos 2 rabbit pAb

Cat No.: ES7253

For research use only

Overview

Product Name Sos 2 rabbit pAb

Host species Rabbit

Applications WB;ELISA;IHC

Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000

Immunogen The antiserum was produced against synthesized

peptide derived from human SOS2. AA range:631-680

Specificity Sos 2 Polyclonal Antibody detects endogenous levels

of Sos 2 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Son of sevenless homolog 2

Gene Name SOS2

Cellular localization intracellular, cytosol,

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 150kD
Human Gene ID 6655
Human Swiss-Prot Number Q07890

Alternative Names SOS2; Son of sevenless homolog 2; SOS-2

Background This gene encodes a regulatory protein that is

involved in the positive regulation of ras proteins. Mutations in this gene are associated with Noonan


Syndrome-9. [provided by RefSeq, Jul 2016],




+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com






Western blot analysis of SOS2 Antibody. The lane on the right is blocked with the SOS2 peptide.



Western blot analysis of the lysates from HT-29 cells using SOS2 antibody.



+86-27-59760950

Immunohistochemical analysis of paraffin-embedded human oophoroma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

