

BUB1 rabbit pAb

Cat No.:ES7356

For research use only

Overview

Product Name	BUB1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human BUB1. AA
	range:781-830
Specificity	BUB1 Polyclonal Antibody detects endogenous
	levels of BUB1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Mitotic checkpoint serine/threonine-protein kinase
	BUB1
Gene Name	BUB1
Cellular localization	Nucleus. Chromosome, centromere, kinetochore.
	Nuclear in interphase cells. Accumulates gradually
	during G1 and S phase of the cell cycle, peaks at
	G2/M, and drops dramatically after mitosis.
	Localizes to the outer kinetochore. Kinetochore
	localization is
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	120kD
Human Gene ID	699
Human Swiss-Prot Number	O43683
Alternative Names	BUB1; BUB1L; Mitotic checkpoint
	serine/threonine-protein kinase BUB1; hBUB1;



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C




Background

BUB1A

This gene encodes a serine/threonine-protein kinase that play a central role in mitosis. The encoded protein functions in part by phosphorylating members of the mitotic checkpoint complex and activating the spindle checkpoint. This protein also plays a role in inhibiting the activation of the anaphase promoting complex/cyclosome. This protein may also function in the DNA damage response. Mutations in this gene have been associated with aneuploidy and several forms of cancer. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],

Western blot analysis of the lysates from HUVECcells using BUB1 antibody.





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C