

UBF rabbit pAb

Cat No.: ES7473

For research use only

Overview

Product Name UBF rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not

yet tested in other applications.

Immunogen Synthesized peptide derived from UBF . at AA range:

420-500

Specificity UBF Polyclonal Antibody detects endogenous levels

of UBF protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Nucleolar transcription factor 1

Gene Name UBTF

Cellular localization Nucleus, nucleolus.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 89kD
Human Gene ID 7343
Human Swiss-Prot Number P17480

Alternative Names UBTF; UBF; UBF1; Nucleolar transcription factor 1;

Autoantigen NOR-90; Upstream-binding factor 1;

UBF-1

Background This gene encodes a member of the HMG-box

DNA-binding protein family. The encoded protein plays a critical role in ribosomal RNA transcription as a key component of the pre-initiation complex, mediating the recruitment of RNA polymerase I to rDNA promoter regions. The encoded protein may



+86-27-59760950 ELKbio@ELKbiotech.com



also play important roles in chromatin remodeling and pre-rRNA processing, and its activity is regulated by both phosphorylation and acetylation.

Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. Pseudogenes of this gene are located on the short arm of chromosomes 3, 11 and X and the long arm of chromosome 11. [provided by RefSeq, Aug 2011],



+86-27-59760950