

## 14-3-3 ζ/δ (phospho Thr232) rabbit pAb

Cat No.:ES7524

For research use only

## Overview

**Product Name** 14-3-3 ζ/δ (phospho Thr232) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

**Recommended dilutions** Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications. The antiserum was produced against synthesized

The antiserum was produced against synthesized peptide derived from human 14-3-3 zeta/delta around the phosphorylation site of Thr232. AA

range:196-245

**Specificity** Phospho-14-3-3 ζ/δ (T232) Polyclonal Antibody

detects endogenous levels of 14-3-3  $\zeta/\delta$  protein only

when phosphorylated at T232.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Protein Name** 14-3-3 protein zeta/delta

Gene Name YWHAZ

Cellular localization Cytoplasm . Melanosome . Located to stage I to

stage IV melanosomes.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 28kD
Human Gene ID 7534
Human Swiss-Prot Number P63104

Alternative Names YWHAZ; 14-3-3 protein zeta/delta; Protein kinase C

inhibitor protein 1; KCIP-1

**Background** This gene product belongs to the 14-3-3 family of


proteins which mediate signal transduction by




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotec



binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene. [provided by RefSeq, Oct 2008],




Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using 14-3-3 zeta/delta (Phospho-Thr232) Antibody



Immunohistochemistry analysis of paraffin-embedded human brain, using 14-3-3 zeta/delta (Phospho-Thr232) Antibody. The picture on the right is blocked with the phospho peptide.







Western blot analysis of lysates from Jurkat cells treated with UV 15', using 14-3-3 zeta/delta (Phospho-Thr232) Antibody. The lane on the right is blocked with the phospho peptide.



+86-27-59760950