

ADAMTS-12 rabbit pAb

Cat No.:ES7651

For research use only

Overview

Product Name	ADAMTS-12 rabbit pAb
Host species	Rabbit
Applications	IHC;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/20000. Not yet tested in other applications.
Immunogen	Synthesized peptide derived from ADAMTS-12 . at
<u> </u>	AA range: 1100-1180
Specificity	ADAMTS-12 Polyclonal Antibody detects
	endogenous levels of ADAMTS-12 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	A disintegrin and metalloproteinase with
	thrombospondin motifs 12
Gene Name	ADAMTS12
Cellular localization	Secreted, extracellular space, extracellular matrix .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	81792
Human Swiss-Prot Number	P58397
Alternative Names	ADAMTS12; A disintegrin and metalloproteinase
	with thrombospondin motifs 12; ADAM-TS 12;
	ADAM-TS12; ADAMTS-12
Background	This gene encodes a member of the ADAMTS (a
-	disintegrin and metalloproteinase with
	thrombospondin motifs) protein family. Members of
	the family share several distinct protein modules,
	including a propeptide region, a metalloproteinase



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



domain, a disintegrin-like domain, and a thrombospondin type 1 (TS-1) motif. Individual members of this family differ in the number of C-terminal TS-1 motifs, and some have unique C-terminal domains. The enzyme encoded by this gene contains eight TS-1 motifs. It may play roles in pulmonary cells during fetal development or in tumor processes through its proteolytic activity or as a molecule potentially involved in regulation of cell adhesion. [provided by RefSeq, Jul 2008],

Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C