

## SerpinA6 rabbit pAb

## Cat No.:ES7846

For research use only

## Overview

Product Name	SerpinA6 rabbit pAb
Host species	Rabbit
Applications	IHC;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
<b>Recommended dilutions</b>	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/5000. Not yet tested in other applications.
Immunogen	Synthesized peptide derived from the Internal
	region of human SerpinA6.
Specificity	SerpinA6 Polyclonal Antibody detects endogenous
	levels of SerpinA6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Corticosteroid-binding globulin
Gene Name	SERPINA6
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	866
Human Swiss-Prot Number	P08185
Alternative Names	SERPINA6; CBG; Corticosteroid-binding globulin;
	CBG; Serpin A6; Transcortin
Background	This gene encodes an alpha-globulin protein with
	corticosteroid-binding properties. This is the major
	transport protein for glucorticoids and progestins in
	the blood of most vertebrates. The gene localizes to
	a chromosomal region containing several closely
	related serine protease inhibitors which may have
	evolved by duplication events. [provided by RefSeq,



+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



## Jul 2008],



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C