

Myt 1 (phospho Ser83) rabbit pAb

Cat No.:ES7956

For research use only

Overview

Product Name	Myt 1 (phospho Ser83) rabbit pAb
Host species	Rabbit
Applications	IHC;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/5000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human MYT1 around the
	phosphorylation site of Ser83. AA range:49-98
Specificity	Phospho-Myt 1 (S83) Polyclonal Antibody detects
	endogenous levels of Myt 1 protein only when
	phosphorylated at S83.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Membrane-associated tyrosine- and
	threonine-specific cdc2-inhibitory kinase
Gene Name	PKMYT1
Cellular localization	Endoplasmic reticulum membrane ; Peripheral
	membrane protein . Golgi apparatus membrane ;
	Peripheral membrane protein .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9088
Human Swiss-Prot Number	Q99640
Alternative Names	PKMYT1; MYT1; Membrane-associated tyrosine- and
	threonine-specific cdc2-inhibitory kinase; Myt1
	kinase
Background	This gene encodes a member of the
	El Khia@El Khiatash som



+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C




serine/threonine protein kinase family. The encoded protein is a membrane-associated kinase that negatively regulates the G2/M transition of the cell cycle by phosphorylating and inactivating cyclin-dependent kinase 1. The activity of the encoded protein is regulated by polo-like kinase 1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, May 2012],



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using MYT1 (Phospho-Ser83) Antibody

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using MYT1 (Phospho-Ser83) Antibody. The picture on the right is blocked with the phospho peptide.





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C