

DCL-1 rabbit pAb

Cat No.: ES8135

For research use only

Overview

Product Name DCL-1 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human CD302. AA

range:51-100

Specificity DCL-1 Polyclonal Antibody detects endogenous

levels of DCL-1 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name CD302 antigen

Gene Name CD302

Cellular localization Membrane; Single-pass type I membrane protein.

Cell projection, filopodium . Cytoplasm, cell cortex . Cell projection, microvillus . Colocalizes with F-actin in filopodia, cellular cortex and microvilli of the

apical cell surface.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 26kD
Human Gene ID 9936
Human Swiss-Prot Number Q8IX05

Alternative Names CD302; CLEC13A; DCL1; KIAA0022; CD302 antigen;

C-type lectin BIMLEC; C-type lectin domain family 13 member A; DEC205-associated C-type lectin 1; Type



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




Background

I transmembrane C-type lectin receptor DCL-1; CD antigen CD302

CD302 is a C-type lectin receptor involved in cell adhesion and migration, as well as endocytosis and phagocytosis (Kato et al., 2007 [PubMed 17947679]).[supplied by OMIM, Aug 2008],

Immunohistochemistry analysis of paraffin-embedded human lymph node tissue, using CD302 Antibody. The picture on the right is blocked with the synthesized peptide.



-- 117 -- 85 -- 48 -- 34 -- 26 -- 19 (kD) Western blot analysis of lysates from HT-29 cells, using CD302 Antibody. The lane on the right is blocked with the synthesized peptide.

