

Atg9a rabbit pAb

Cat No.:ES8476

For research use only


Overview

Product Name	Atg9a rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000, ELISA 1:10000-20000
Immunogen	Synthetic Peptide of Atg9a
Specificity	Atg9a Polyclonal Antibody detects endogenous levels of Atg9a
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Atg9a
Gene Name	ATG9A APG9L1
Cellular localization	Preautophagosomal structure membrane ; Multi-pass membrane protein . Cytoplasmic vesicle, autophagosome membrane ; Multi-pass membrane protein . Golgi apparatus, trans-Golgi network membrane ; Multi-pass membrane protein . Late endosome membrane ; Multi-p
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	85kD
Human Gene ID	79065
Human Swiss-Prot Number	Q7Z3C6
Alternative Names	ATG9 autophagy related 9 homolog A (S. cerevisiae)
Background	function:Plays a role in autophagy.,similarity:Belongs to the ATG9 family,





ELK Biotechnology



Western Blot analysis of 293T 3T3 cells using Atg9a Polyclonal Antibody diluted at 1:800. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C