

## **GMF-γ** rabbit pAb

Cat No.: ES8725

For research use only

## Overview

Product Name GMF-γ rabbit pAb

Host species Rabbit
Applications IHC;IF;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions IHC-p 1:50-200, ELISA 1:10000-20000

Immunogen Synthetic peptide from human protein at AA range:

1-50

**Specificity** The antibody detects endogenous GMF-γ

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Glia maturation factor gamma (GMF-gamma)

Gene Name GMFG
Cellular localization intracellular,

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

**Human Gene ID** 9535 **Human Swiss-Prot Number** 060234

+86-27-59760950

Alternative Names Glia maturation factor gamma (GMF-gamma)


**Background** similarity:Belongs to the actin-binding proteins ADF

family. GMF subfamily., similarity: Contains 1 ADF-H domain., tissue specificity: Expressed predominantly


in lung, heart and placenta.,







Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

