

AHNK rabbit pAb

Cat No.:ES8947

For research use only

Overview

Product Name	AHNK rabbit pAb
Host species	Rabbit
Applications	IHC;IF
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	IHC-p 1:50-300
Immunogen	Synthesized peptide derived from human protein .
	at AA range: 5720-5800
Specificity	AHNK Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Neuroblast differentiation-associated protein
	AHNAK (Desmoyokin)
Gene Name	AHNAK PM227
Cellular localization	Nucleus.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	647kD
Human Gene ID	79026
Human Swiss-Prot Number	Q09666
Alternative Names	
Background	function:May be required for neuronal cell
	differentiation., similarity: Contains 1 PDZ (DHR)
	domain.,subunit:Interacts with DYSF; the interaction
	is direct and Ca(2+)-independent.,




ELKbio@ELKbiotech.com

h.com ww

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Tris-EDTA,pH9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200(4° overnight.3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C