

## **PGFRL** rabbit pAb

## Cat No.:ES9063

For research use only

## Overview

Product Name	PGFRL rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
<b>Recommended dilutions</b>	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein .
	at AA range: 50-130
Specificity	PGFRL Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Platelet-derived growth factor receptor-like protein
	(PDGFR-like protein) (PDGF receptor beta-like tumor
	suppressor)
Gene Name	PDGFRL PRLTS
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	41kD
Human Gene ID	5157
Human Swiss-Prot Number	Q15198
Alternative Names	
Background	This gene encodes a protein with significant
	sequence similarity to the ligand binding domain of
	platelet-derived growth factor receptor beta.
	Mutations in this gene, or deletion of a
	chromosomal segment containing this gene, are
	associated with sporadic hepatocellular carcinomas,
	colorectal cancers, and non-small cell lung cancers.




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





This suggests this gene product may function as a tumor suppressor. [provided by RefSeq, Jul 2008],

Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C