

SMG9 rabbit pAb

Cat No.: ES9092

For research use only

Overview

Product Name SMG9 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 70-150

Specificity SMG9 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Protein SMG9 (Protein smg-9 homolog)

Gene Name SMG9 C19orf61 **Cellular localization** intracellular,cytosol,

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 57kD
Human Gene ID 56006
Human Swiss-Prot Number Q9H0W8

Alternative Names

Background SMG9, nonsense mediated mRNA decay

factor(SMG9) Homo sapiens This gene encodes a regulatory subunit of the SMG1 complex, which plays a critical role in nonsense-mediated mRNA decay (NMD). Binding of the encoded protein to the SMG1 complex kinase scaffold protein results in the inhibition of its kinase activity. Mutations in this gene cause a multiple congenital anomaly syndrome

in human patients, characterized by brain




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com





malformation, congenital heart disease and other features. [provided by RefSeq, Jul 2016],

Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000



+86-27-59760950