

## **ALKB7** rabbit pAb

Cat No.: ES9269

For research use only

## Overview

Product Name ALKB7 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 10-90

Specificity ALKB7 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Probable alpha-ketoglutarate-dependent

dioxygenase ABH7 (EC 1.14.11.-) (Alkylated DNA repair protein alkB homolog 7) (Spermatogenesis

cell proliferation-related protein) (Spermatogenesis-associated pr

Gene Name ALKBH7 ABH7 SPATA11 UNQ6002/PRO34564

**Cellular localization** Mitochondrion matrix.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.


Clonality Polyclonal
Concentration 1 mg/ml
Observed band 24kD
Human Gene ID 84266
Human Swiss-Prot Number Q9BT30

**Alternative Names** 

**Background** similarity: Belongs to the alkB family.,







Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

