

ACHG rabbit pAb

Cat No.: ES9337

For research use only

Overview

Product Name ACHG rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 30-110

Specificity ACHG Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Acetylcholine receptor subunit gamma

Gene Name CHRNG ACHRG

Cellular localization Cell junction, synapse, postsynaptic cell membrane;

Multi-pass membrane protein. Cell membrane;

Multi-pass membrane protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 56kD
Human Gene ID 1146
Human Swiss-Prot Number P07510

Alternative Names

Background The mammalian muscle-type acetylcholine receptor

is a transmembrane pentameric glycoprotein with two alpha subunits, one beta, one delta, and one epsilon (in adult skeletal muscle) or gamma (in fetal and denervated muscle) subunit. This gene, which encodes the gamma subunit, is expressed prior to the thirty-third week of gestation in humans. The



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



gamma subunit of the acetylcholine receptor plays a role in neuromuscular organogenesis and ligand binding and disruption of gamma subunit expression prevents the correct localization of the receptor in cell membranes. Mutations in this gene cause Escobar syndrome and a lethal form of multiple pterygium syndrome. Muscle-type acetylcholine receptor is the major antigen in the autoimmune disease myasthenia gravis.[provided by RefSeq, Sep 2009],



+86-27-59760950