

## B3A2 rabbit pAb

Cat No.: ES9393

For research use only

## Overview

Product Name B3A2 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 830-910

**Specificity** B3A2 Polyclonal Antibody detects endogenous levels

of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Anion exchange protein 2 (AE 2) (Anion exchanger 2)

(Non-arthreid band 3 like protein) (RND31) (Solute)

(Non-erythroid band 3-like protein) (BND3L) (Solute

carrier family 4 member 2)

Gene Name SLC4A2 AE2 EPB3L1 HKB3 MPB3L

Cellular localizationMembrane; Multi-pass membrane protein.PurificationThe antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 136kD
Human Gene ID 6522
Human Swiss-Prot Number P04920

**Alternative Names** 

+86-27-59760950

Background alternative products: Additional isoforms seem to

exist, function: Plasma membrane anion exchange protein of wide distribution., similarity: Belongs to

the anion exchanger (TC 2.A.31) family.,

