

## ABCD3 rabbit pAb

Cat No.: ES9425

For research use only

## Overview

Product Name ABCD3 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 40-120

Specificity ABCD3 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name ATP-binding cassette sub-family D member 3 (70 kDa

peroxisomal membrane protein) (PMP70)

Gene Name ABCD3 PMP70 PXMP1

Cellular localization Peroxisome membrane ; Multi-pass membrane

protein.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 72kD
Human Gene ID 5825
Human Swiss-Prot Number P28288

**Alternative Names** 

**Background** The protein encoded by this gene is a member of

the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



the ALD subfamily, which is involved in peroxisomal import of fatty acids and/or fatty acyl-CoAs in the organelle. All known peroxisomal ABC transporters are half transporters which require a partner half transporter molecule to form a functional homodimeric or heterodimeric transporter. This peroxisomal membrane protein likely plays an important role in peroxisome biogenesis. Mutations have been associated with some forms of Zellweger syndrome, a heterogeneous group of peroxisome assembly disorders. Alternative splicing results in multiple transcript variants



+86-27-59760950