

ABCF3 rabbit pAb

Cat No.: ES9426

For research use only

Overview

Product Name ABCF3 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 40-120

Specificity ABCF3 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name ATP-binding cassette sub-family F member 3

Gene Name ABCF3

Cellular localization cell-cell adherens junction, membrane,

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 77kD
Human Gene ID 55324
Human Swiss-Prot Number Q9NUQ8

Alternative Names

+86-27-59760950

Background caution:Lacks transmembrane domains and is

probably not involved in

transport., similarity: Belongs to the ABC transporter family. EF3 subfamily., similarity: Contains 2 ABC

transporter domains.,

