

## CTNA2 rabbit pAb

## Cat No.:ES9508

For research use only

## Overview

Product Name	CTNA2 rabbit pAb	
Host species	Rabbit	
Applications	WB;ELISA	
Species Cross-Reactivity	Human; Mouse	
<b>Recommended dilutions</b>	WB 1:500-2000 ELISA 1:5000-20000	
Immunogen	Synthesized peptide derived from part region of	
	human protein	
Specificity	CTNA2 Polyclonal Antibody detects endogenous	
	levels of protein.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20 $^\circ\!\mathrm{C}$ . Avoid repeated freeze-thaw cycles.	
Protein Name	Catenin alpha-2 (Alpha N-catenin)	
	(Alpha-catenin-related protein)	
Gene Name	CTNNA2 CAPR	
<b>Cellular localization</b>	Cell membrane ; Peripheral membrane protein ;	
	Cytoplasmic side . Cytoplasm . Cytoplasm,	
	cytoskeleton . Cell junction, adherens junction . Cell	
	projection, axon . Nucleus .	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band	104kD	
Human Gene ID	1496	
Human Swiss-Prot Number	P26232	
Alternative Names		
Background	similarity:Belongs to the vinculin/alpha-catenin	
	family.,	
		<i></i> ////////////////////////////////



ELKbio@ELKbiotech.com

om www.elk

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C