

PDE3B rabbit pAb

Cat No.:ES9529

For research use only

Overview

Product Name	PDE3B rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of
-	human protein
Specificity	PDE3B Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
Protein Name	cGMP-inhibited 3',5'-cyclic phosphodiesterase B (EC
	3.1.4.17) (CGIPDE1) (CGIP1) (Cyclic GMP-inhibited
	phosphodiesterase B) (CGI-PDE B)
Gene Name	PDE3B
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	122kD
Human Gene ID	5140
Human Swiss-Prot Number	Q13370
Alternative Names	
Background	catalytic activity:Nucleoside 3',5'-cyclic phosphate +
	H(2)O = nucleoside 5'-phosphate.,enzyme
	regulation:Inhibited by cGMP.,function:May play a
	role in fat metabolism., similarity: Belongs to the
	cyclic nucleotide phosphodiesterase family.,tissue
	specificity:Abundant in adipose tissues.,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elk

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C