

CNTFR rabbit pAb

Cat No.:ES9543

For research use only

Overview

Product Name	CNTFR rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human; Mouse; Rat
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of
	human protein
Specificity	CNTFR Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Ciliary neurotrophic factor receptor subunit alpha
	(CNTF receptor subunit alpha) (CNTFR-alpha)
Gene Name	CNTFR
Cellular localization	Cell membrane; Lipid-anchor, GPI-anchor.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	40kD
Human Gene ID	1271
Human Swiss-Prot Number	P26992
Alternative Names	
Background	This gene encodes a member of the type 1 cytokine
	receptor family. The encoded protein is the
	ligand-specific component of a tripartite receptor
	for ciliary neurotrophic factor, which plays a critical
	role in neuronal cell survival, differentiation and
	gene expression. Binding of ciliary neurotrophic
	factor to the encoded protein recruits the
	transmembrane components of the receptor, gp130



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



and leukemia inhibitory factor receptor, facilitating signal transduction. Single nucleotide polymorphisms in this gene may be associated with variations in muscle strength, as well as early onset of eating disorders. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, May 2011],

Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C