

CCDC8 rabbit pAb

Cat No.:ES9557

For research use only

Overview

Product Name CCDC8 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity CCDC8 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Coiled-coil domain-containing protein 8

Gene Name CCDC8

Cellular localization Cytoplasm . Cytoplasm, cytoskeleton, microtubule

organizing center, centrosome.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 59kD
Human Gene ID 83987
Human Swiss-Prot Number Q9H0W5

Alternative Names

Background This gene encodes a coiled-coil domain-containing

protein. The encoded protein functions as a cofactor required for p53-mediated apoptosis following DNA damage, and may also play a role in growth through interactions with the cytoskeletal adaptor protein obscurin-like 1. Mutations in this gene are a cause of 3M syndrome-3 (3M3). [provided by RefSeq, Dec

2011],

+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

