

## FBLN1 rabbit pAb

Cat No.: ES9659

For research use only

## Overview

Product Name FBLN1 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

**Specificity** FBLN1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at  $-20^{\circ}$ C. Avoid repeated freeze-thaw cycles.

Protein Name Fibulin-1 (FIBL-1)
Gene Name FBLN1 PP213

Cellular localizationSecreted, extracellular space, extracellular matrix.PurificationThe antibody was affinity-purified from rabbit<br/>antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 77kD
Human Gene ID 2192
Human Swiss-Prot Number P23142

**Alternative Names** 

**Background** Fibulin 1 is a secreted glycoprotein that becomes

incorporated into a fibrillar extracellular matrix. Calcium-binding is apparently required to mediate its binding to laminin and nidogen. It mediates platelet adhesion via binding fibrinogen. Four splice

variants which differ in the 3' end have been

identified. Each variant encodes a different isoform, but no functional distinctions have been identified among the four variants. [provided by RefSeq, Jul



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



2008],



+86-27-59760950