

## **QSER1** rabbit pAb

Cat No.: ES9675

For research use only

## Overview

Product Name QSER1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

**Species Cross-Reactivity** Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

**Specificity** QSER1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at  $-20^{\circ}$ C. Avoid repeated freeze-thaw cycles.

Protein Name Glutamine and serine-rich protein 1

Gene Name QSER1

Cellular localization Chromosome .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 190kD
Human Gene ID 79832
Human Swiss-Prot Number Q2KHR3

**Alternative Names** 

+86-27-59760950

Background PTM:Phosphorylated upon DNA damage, probably

by ATM or ATR.,



ELKbio@ELKbiotech.com