

## NDE1 rabbit pAb

Cat No.: ES9931

For research use only

## Overview

Product Name NDE1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Rat; Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 270-350

**Specificity** NDE1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Nuclear distribution protein nudE homolog 1 (NudE)

Gene Name NDE1 NUDE

**Cellular localization** Cytoplasm, cytoskeleton. Cytoplasm, cytoskeleton,

microtubule organizing center, centrosome.

Chromosome, centromere, kinetochore. Cytoplasm, cytoskeleton, spindle. Cleavage furrow. Localizes to the interphase and S phase centrosome. During

mitosis, partia

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 38kD
Human Gene ID 54820
Human Swiss-Prot Number Q9NXR1

**Alternative Names** 


**Background** This gene encodes a member of the nuclear

distribution E (NudE) family of proteins. The

encoded protein is localized at the centrosome and interacts with other centrosome components as







part of a multiprotein complex that regulates dynein function. This protein plays an essential role in microtubule organization, mitosis and neuronal migration. Mutations in this gene cause lissencephaly 4, a disorder characterized by lissencephaly, severe brain atrophy, microcephaly, and severe mental retardation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2012],

Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C