

## P210L rabbit pAb

Cat No.: ES9943

For research use only

## Overview

Product Name P210L rabbit pAb

Host species Rabbit
Applications WB;ELISA

**Species Cross-Reactivity** Human;Rat;Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 980-1060

**Specificity** P210L Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20 °C. Avoid repeated freeze-thaw cycles.

Protein Name Nuclear pore membrane glycoprotein 210-like

(Nucleoporin 210 kDa-like) (Nucleoporin

Nup210-like)

Gene Name NUP210L

Cellular localizationMembrane ; Single-pass membrane protein .PurificationThe antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band207kDHuman Gene ID91181Human Swiss-Prot NumberQ5VU65

**Alternative Names** 

**Background** similarity: Belongs to the NUP210 family.,

