



ELK Biotechnology
For research use only.

dNTP Mix (PCR Grade) 10 mM each

Catalog No.	Specification	Storage/Shelf life
EQ011-1	1 mL	-20°C/1 year
EQ011-2	0.1 mL	-20°C/1 year

Introductoin

This product has undergone rigorous quality testing with a purity of over 99% and does not contain RNase and DNase. dNTPs Mix is an equimolar mixture of dATP, dGTP, dCTP and dTTP. When used as a substrate for PCR reaction, it can be diluted to a suitable concentration with sterile ultrapure water.