

Recombinant Human LIF

Catalog #	EPT153
Expression Host	E.coli
DESCRIPTION	Recombinant Human Leukemia Inhibitory Factor is
	produced by our E.coli expression system and the
	target gene encoding Ser23-Phe202 is expressed.
Accession	P15018
Synonyms	Leukemia Inhibitory Factor; LIF;
	Differentiation-Stimulating Factor; D Factor;
	Melanoma-Derived LPL Inhibitor; MLPLI; Emfilermin;
	LIF; HILDA
Mol Mass	19.7 KDa
AP Mol Mass	18 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.001 ng/ μ g (0.01 EU/ μ g) as determined by
	LAL test.
FORMULATION	Lyophilized from a 0.2 μ m filtered solution of PBS, pH
	7.4.



+86-27-59760950

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C

www.elkbiotech.com

ELKbio@ELKbiotech.com



RECONSTITUTION Always centrifuge tubes before opening.Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. SHIPPING The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below. Lyophilized protein should be stored at < -20 °C, STORAGE though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at $< -20^{\circ}$ C for 3 months. BACKGROUND Leukemia Inhibitory Factor (LIF) is a lymphoid factor that promotes long-term maintenance of embryonic stem cells by suppressing spontaneous differentiation. LIF has a number of other activities including cholinergic neuron differentiation, control of stem cell pluripotency, bone and fat metabolism, mitogenesis



+86-27-59760950

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.G

www.elkbiotech.com

ELKbio@ELKbiotech.com



of certain factor dependent cell lines and promotion of megakaryocyte production in vivo. Human and murine mature LIF exhibit a 78% sequence identity at the amino acid level. Human LIF is equally active on human and mouse cells. Murine LIF is approximately 1000 fold less active on human cells than human LIF.



SDS-PAGE



+86-27-59760950 ELKbio@ELKbiotech.com

otech.com www

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C