

Recombinant Human CD32a

(H167,C-6His)

Catalog #	EPT162
-----------	--------

Expression Host Human Cells

DESCRIPTION Recombinant Human Low Affinity Immunoglobulin Gamma Fc Region Receptor II-A(H167) is produced by our Mammalian expression system and the target gene encoding Ala36-IIe218 is expressed with a 6His tag at the C-terminus. It is identical to FCGR2A131H/R in the reference.

Accession P12318

Synonyms Low affinity immunoglobulin gamma Fc region receptor II-a; IgG Fc receptor II-a; CDw32; Fc-gamma RII-a; Fc-gamma-RIIa; FcRII-a; CD32; FCGR2A; FCG2; FCGR2A1; IGFR2

Mol Mass 21.1 KDa

AP Mol Mass 25-32 KDa, reducing conditions

 Purity
 Greater than 95% as determined by reducing

 SDS-PAGE.



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com 23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Less than 0.1 ng/ μ g (1 EU/ μ g) as determined by LAL Endotoxin test. FORMULATION Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. RECONSTITUTION Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. SHIPPING The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below. STORAGE Lyophilized protein should be stored at < -20 °C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20° C for 3 months. BACKGROUND Human Fc y Rs are divided into three classes designated FcyRI (CD64), FcyRII (CD32), and FcyRIII



+86-27-59760950

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.G

www.elkbiotech.com

ELKbio@ELKbiotech.com



(CD16), which generate multiple isoforms, are recognized. The activating - type receptor either has or associates non-covalently with an accessory subunit that has an immunoreceptor tyrosine-based activation motif (ITAM) in its cytoplasmic domain. FcyRI binds IgG with high affinity and functions during early immune responses, whereas FcyRII and RIII are low affinity receptors that recognize IgG as aggregates surrounding multivalent antigens during late immune responses.Human CD32, also known as Low affinity immunoglobulin y Fc region receptor II-a (IgG Fc receptor II-a), FcyRII A or FCGR2A Protein, is expressed on cells of both myeloid and lymphoid lineages as well as on cells of non-hematopoietic origin. Associated with an ITAM-bearing adapter subunit, FcRy, CD32a (Fc y RII A) delivers an activating signal upon ligand binding, and results in the initiation of inflammatory including responses cytolysis, phagocytosis, degranulation, and cytokine production. The responses can be modulated by signals from the co-expressed inhibitory receptors such as Fcy RII B, and the strength of the signal is dependent on the



+86-27-59760950 ELK

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



ratio of expression of the activating and inhibitory receptors.



SDS-PAGE



+86-27-59760950 ELKbio

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C